


Calibration options

Coordinate axes
Calibration stick/tape
Offset origin
Calibration points
Circle fitter

Origin	Angle	Scale
X	X	
	X	X
X		
X	X	X
X		



Coordinate axes




- Most common way to set origin and angle
- Drag origin or nudge with arrow keys
- Set origin pixel position on toolbar
- Drag x-axis to tilt, shift-drag for 5° steps
- Set x-axis angle on toolbar

Tracker Help: coordinate axes



Calibration stick and tape




- Most common way to set scale
- Drag either end or nudge with arrow keys
- Dbl-click length readout to set new value
- Set length value on toolbar
- Set x-axis angle on toolbar
- · Stick and tape behave differently when dragging end

Tracker Help: calibration stick and tape



Offset origin or single calibration point




- · Handy when origin is off screen
- Drag or nudge to move origin
- Set coordinates on toolbar to move origin

Tracker Help: offset origin



Calibration points (pair)




- Most powerful calibration tool
- To use, position both points at known locations
- · Set coordinates on toolbar to set origin, angle, tilt
- Drag or nudge either point to adjust

Tracker Help: calibration points

DL browser: fluor\_lamp\_3000K.zip



Circle fitter



- Fit circle to 3+ points (may be point mass steps)
- · Choose "move origin to center" from popup or track menu
- Useful for circular motion, esp when origin off screen

Tracker Help: circle fitter

DL browser: Rotating fan.mp4